

THE MULTI-PURPOSE BANDAGE



Trauma/Wound Dressing

The Multi-Purpose Bandage is the **first major innovation** of its kind in over 50 years!

First Care Products Ltd. has developed a new internationally patented, first-aid device for the staunching of blood flow from traumatic hemorrhagic wounds in pre-hospital emergency situations.

The Multi-Purpose Bandage is an ideal solution for emergency treatment, where every second counts. *The product consolidates many first-aid devices such as: primary dressing, pressure applicators, secondary dressing and a closure bar into one unit.*

The Multi-Purpose Bandage is easy for the professional and non-professional care-giver to apply and perform multiple functions.

The Multi-Purpose Bandage was developed with the most extreme conditions in mind. Designed to treat every possible bleeding wound, the sterile, non-adherent bandage applies pressure to any site, is easily wrapped, secured and acts as a tourniquet in cases of severe bleeding. *The product is so easy to use that even an injured person can apply the bandage with only one hand.*



APPLICATIONS:

- Home first-aid kit
- Automotive first-aid kit
- Commercial and Industrial sites
- First Responders
- Emergency Medical Services
- Military

THE MULTI-PURPOSE BANDAGE: Advantages

- Immediate direct pressure
- Easy personal and professional application
- Consolidation into a single unit:
 - **primary dressing**
 - **pressure applicators**
 - **secondary dressing**
 - **closure bar**
- Significant time and cost savings

FIRST CARE PRODUCTS

3050 Post Oak Blvd., Suite 1710 • Houston, TX 77056 • Tel.: 713-723-6000 • Fax: 713-723-6221

email: info@firstcareproducts.com

Website: www.firstcareproducts.com